



Armschoner BioClean-C

pure¹¹-Nr.: 03504, Hersteller: BioClean





Zusammenfassung

- Neue pure11-Artikelnummer (ab 01.07.2023): 1103504
- Material: Clean-Tough
- Gammabestrahlt
- Mit elastischem Gummibund am Oberarm und Polyester-Trikotmanschette am unteren Ende
- Für bequemen Sitz
- Maximaler Komfort und Schutz
- Zyto-geeignet
- Einzeln verpackt

Empfohlene Reinraumklassen

ISO 3 4 5 6 7 8 9

GMP A/B C D

Produktvarianten

pure¹¹-Nr.: 03504B

Farbe: Blau / Herst.-Nr.: S-BCSC / VE: 90 Paar

Quelle: https://www.pure11.de/armschoner-bioclean-c-0



BIOCLEAN-C CHEMOTHERAPY Sterile Protective Sleeve Covers S-BCSC

About

The BioClean-C Sterile Chemotherapy
Protective Sleeve Covers have been
specially developed for protection against
a range of chemotherapy drugs** and to
be compatible with ISO Class 4
cleanrooms. Constructed from
lightweight low-linting CleanTough blue
material, the sleeve covers feature
elasticated arms for stability on the arm
and has been tailored for quick and
simple donning.

Specifications

MATERIAL:

BioClean-C CleanTough blue material & 100% polyester cuffs

COLOUR:

Blue

CONSTRUCTION:

Ultrasonically bonded seams covered with protective tape

CLASSIFICATION:

PPE Cat 3 Type PB [6]

CHARACTERISTICS: Low particulating

COMPATIBILITY:

ISO Class 4

STERILIZATION:
Gamma irradiation

Features

- Ultrasonically bonded seams with protective tape
- Elasticated for secure fit
- Lightweight low-linting CleanTough material
- Tested against permeation standard ASTM F739-12
- Silicone-free







Quality Standards

- Manufactured in a facility holding ISO 9001:2015 and (EU)2016/425
 Module D certifications
- Conforms to Category 3 Personal Protective Equipment Regulation (EU) 2016/425
- Tested against standard EN ISO 13034:2005 + A1:2009 protective clothing against liquid chemicals (Type PB [6])
- Tested against standard EN 14325 protective clothing against chemicals
- Tested against permeation standard ASTM F739-12
- Processed in an NEBB certified ISO Class 4 environment
- Tested against standards: EN 530, Abrasion resistance; EN 6530, Resistance to penetration by liquids; EN 13934-1, Tensile strength; EN 9073-4, Trapezoidal tear resistance; EN 863, Puncture resistance; and EN 13935-2, Seam strength.
- Sterilized in accordance with ISO 11137-1:2006

Sterilization

Method: Gamma irradiation Minimum Dose: 25kGy Sterility Assurance Level: 10⁻⁶

Shelf Life & Storage

Three (3) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

Performance Test Results - BioClean-C CleanTough Blue Material

TEST	RESUL	г	PERFORMANCE CLASS (EN 14325)	
Abrasion Resistance	10 to 100 d	ycles	1	
Trapezoidal Tear Resistance	CD 32.3	N	2	
Trapezoidal Tear Resistance	MD 50.6	N	3	
Tensile Strength	CD 47 N		1	
Tensile Strength	MD 94	N	2	
Puncture Resistance	10 N		1	
Repellence to Liquids	30% H ₂ SO ₄	99.0%	3	
	10% NaOH	98.4%	3	
	O-Xylene	91.4%	1	
	Butan-1-ol	93.9%	3	
Penetration by Liquids	30% H ₂ SO ₄	0%	3	
	10% NaOH	0%	3	
	O-Xylene	0%	3	
	Butan-1-ol	0%	3	
Seam Strength ¹	74 N		2	

Particle Shedding Test Results

TEST	RESULT	CATEGORY (IEST-RP-CC003.3)	
Particle Shedding (Helmke Drum Test)	<260	Category I	





BioClean-C™ CleanTough Blue Material - ASTM F739-12 Test Method Results

Mean Breakthrough Time (MBT), Minutes

Breakthrough of the test chemical is deemed to have occurred when the permeation rate has reached 0.1 μ g/cm²/min

CISPLATIN	CARMUSTINE	CYCLO- PHOSPHA- MIDE	DOXO- RUBICIN HYDRO- CHLORIDE	5- FLUORO- URACIL	METHO- TREXATE	ETOPOSIDE	PACLITAXEL	THIOTEPA
>240	>240	>240	>240	>240	>240	<6	>240	55 (46,88,31)

Results achieved under controlled laboratory conditions, by accredited external testing laboratory.

TO ORDER

Re-OrderChestHeightSizeS-BCSCN/AN/AUniversal

Packing

One pair per sealed inner PE bag; 15 inner bags per sealed outer PE bag; six outer bags per lined carton (90 pairs)

For additional information visit us at www.bioclean.com, or call us at +44 1638 663338

Nitritex Ltd

Minton Enterprise Park Oaks Drive Newmarket, Suffolk CB8 7YY United Kingdom

S-BCSC/1905/PDS33



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of aprons/sleeve covers for a specific application.

**Please see product validation pack or contact Ansell customer service for specific data on use of garments with cytotoxic drugs. Garments used for protection against such drugs must be selected specifically for the type of chemicals used.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2019 Ansell Limited. All Rights Reserved.