CELOS[®]

a pure¹¹ brand

CELOS GLOVE 1.10

English
German
Dutch
Spanish
Portuguese
French
Danish
Polish
Czech
Hungarian
Macedonian
Slovakian

Instruction manual

EN Nitrile Glove CELOS Glove 1.10

Nitrile Glove for Cleanroom Applications

NOTE: THIS PACKAGE INSERT SHOULD BE PROVIDED TO USERS OF THESE GLOVES FOR SAFETY REASONS.

This product is certified as a Category III PPE product in accordance with regulation (EU) 2016/425 EEC.

TÜV Rheinland LGA Products GmbH · Tillystraße 2 · D-90431 Nürnberg · Identification number 0197

Nitrile · White · Ambidextrous · Powder-free · Latex-free · Acceptable Quality Level (AQL): 1.5 · Textured fingertips · Accelerator-free · Low particulate and ion residue levels · Washed with 0.2 µm filtered DI water · Double packaging

USE: These gloves are designed to protect the user's hands from low chemical risks and microbial contamination. Gloves are part of protective clothing and protect the user from potentially toxic substances or organisms. On the other hand, they also serve to protect products and samples from contamination of any kind. These gloves provide no protection against mechanical risks. The testing and labeling of the gloves comply with the PPE Regulation 2016/425, which replaces 89/686/EEC, and the applicable harmonized European standards. Please ensure that the gloves are used exclusively for their intended applications.

NOTE: This information does not indicate the actual duration of protection in the workplace or the differentiation between mixed and pure chemicals. Resistance to chemicals has been assessed under laboratory conditions on samples taken from the inside of the hand only and refers only to the chemicals tested. Verification that the gloves are suitable for the intended use is recommended, as workplace conditions may be different from those of the type test depending on temperature, abrasion, and degradation. If protective gloves have already been used, they may offer less resistance to hazardous chemicals due to changes in their physical attributes. Degradation caused by contact with chemicals, movements, friction, can significantly reduce the actual application time. With aggressive chemicals, degradation may be the most important factor to consider when selecting chemical resistant gloves. According to standard EN ISO 374-5: The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen.

STORAGE RECOMMENDATION: Do not store cleanroom nitrile gloves in locations where the temperature may exceed 104°F (40°C). Opened cleanroom nitrile gloves must be protected from direct sunlight or prolonged fluorescent lighting to prevent discoloration. Incorrect storage of cleanroom nitrile gloves will result in a reduced shelf life and will adversely affect their properties.

ATTENTION: Primary material nitrile gloves do not contain natural rubber latex, ensuring a low allergy risk. Do not use in case of intolerance to nitrile. Inspect gloves for damage before wearing. If there are holes, tears or re-wearing, the protective performance can no longer be guaranteed. Do not use if damaged.

DISPOSAL INFORMATION: The product must be disposed of in accordance with the rules and regulations in force in your country. Gloves which have been contaminated by chemical substances must be disposed of in accordance with the regulations for the chemicals concerned.

You can read the conformity assessment and the issuance of the EU declaration of conformity here: <https://www.pure11.de/de/doc>

Importer: pure11 GmbH · Bavariafilmpatz 7 · 82031 Gruenwald · Germany

No	Chemical Testing	Cas Number	Chemical Code
1	N-heptane	142-82-5	J
2	Hydrogen peroxide	7722-84-1	P
3	Formaldehyde	50-00-0	T
4	Sodium hydroxide	1310-73-2	K
5	Ethanol	64-17-5	-/-
6	Hydrochloric acid	7647-01-0	-/-
7	Sulfuric acid	7664-93-9	-/-

Permeation Test EN16523-1:2015			Degradation Test EN 374-4:2014
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
N-heptane	> 30	Level 2	-40.7

Permeation Test EN16523-1:2015			Degradation Test EN 374-4:2014
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Hydrogen peroxide	> 480	Level 6	31.6

Permeation Test EN16523-1:2015			Degradation Test EN 374-4:2014
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Formaldehyde	> 30	Level 2	18.8

Permeation Test EN16523-1:2015			Degradation Test EN 374-4:2014
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Sodium hydroxide	> 480	Level 6	-24.3

Permeation Test EN16523-1:2015			Degradation Test EN 374-4:2014
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Ethanol (70 %)	> 10	Level 1	n.b.

Permeation Test EN16523-1:2015			Degradation Test EN 374-4:2014
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Hydrochloric acid (30 %)	> 480	Level 6	n.b.

Permeation Test EN16523-1:2015			Degradation Test EN 374-4:2014
Chemical	Breakthrough Time (min.)	Performance Level	Degradation %
Sulfuric acid (50 %)	> 480	Level 6	n.b.

Please consider: This information does not necessarily reflect the actual breakthrough time in the workplace. Chemical resistances have been assessed under laboratory conditions.

Product Code	Size
1105050XS	EX-SMALL Hand Size (5.5)
1105050S	SMALL Hand Size (6.0-6.5)
1105050M	MEDIUM Hand Size (7.0-7.5)
1105050L	LARGE Hand Size (8.0-8.5)
1105050XL	EX-LARGE Hand Size (9.0-9.5)
1105050XXL	EX-EX-LARGE Hand Size (10.0-10.5)
11050503XL	EX-EX-EX-LARGE Hand Size (11.0-11.5)

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Importér: pure11 GmbH · Bavariafilmplatz 7 · 82031 Grünwald · Nemecko