Contec Sterile HydroPure

pure¹¹-Nr.: 09032, Hersteller: CONTEC

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Contec HydroPure

Stabilised 6% hydrogen peroxide

Contec HydroPure is a blend of 6% hydrogen peroxide and water for injection or purified water. Efficacious against bacteria, fungi, moulds, yeasts and spores, Contec HydroPure leaves little to no residue and is available with a guaranteed endotoxin level of less than 0.25 EU/ml making it ideal for use in product contact areas.

Contec HydroPure is 0.2 micron filtered, and filled under unidirectional air flow in a Grade B cleanroom. The sterile product is filled into pre-irradiated components.

Available as either a 1L trigger spray for small areas or in 0.5L, 1L and 5L capped containers for mop head saturation or when larger volumes are required.

Contec Sterile HydroPure 1L trigger spray is blended with water for injection.

Provided either double or triple bagged, the product is designed for ease of entry into pharmaceutical cleanrooms. The 1L trigger sprays are fitted with a protective system to protect the contents during use.

FEATURES	BENEFITS
Hydrogen peroxide breaks down to water and oxygen	• Very low residue so suitable for product contact areas
Guaranteed endotoxin level of less than 0.25 EU/ml	• Safe to use in product contact areas
Filtered to 0.2 microns under Grade A air-flow in a Grade B cleanroom	• Ensures the product is free from contamination and particulates
Sporicidal in 15 mins	• Sufficient sporicidal activity for a Grade A / B area
Sterile version available	• Suitable for use in Grade A and B cleanrooms
Trigger spray can be set to jet or spray	• Large droplet size reduces the risk of inhalation and provides good surface coverage
Triple / double bagged in linear tear packaging	• Each bag is easy to open even when wearing gloves • Facilitates transfer disinfection into isolators and RABS
Not classed as corrosive	• Can be used safely in all areas of the cleanroom

For more information or to request a sample,
e-mail infoeu@contecinc.com or phone + 33 (0) 2 97 43 76 98.

Copyright© 2020 Contec, Inc. All rights reserved. INT109 280220

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

Sterile HydroPure Specification

- Trigger spray - 6% hydrogen peroxide in water for injection
- Capped - 6% hydrogen peroxide in purified water
- 1L trigger spray and 5L capped - less than 0.25 EU/ml
- Sterile filled into pre-irradiated components
- Trigger spray - triple bagged in polyethylene linear tear bags
- Capped - double bagged in polyethylene linear tear bags
- RTU with a 2 year shelf life

Filtered HydroPure Specification

- 6% hydrogen peroxide in purified water
- Filtered to 0.2 microns in a Grade B cleanroom
- Double bagged in polyethylene linear tear bags
- RTU with a 2 year shelf life

Certificates of analysis provide with every batch and full supporting documentation is available on request.

Efficacy

Test	Description	Log Reduction	Time	Test	Description	Log Reduction	Time
EN1276	<i>E. hirae</i>	>5.3	5 min	EN1650	<i>A. brasiliensis</i>	>4.0	15 min
EN1276	<i>S. aureus</i>	>5.4	5 min	EN13704	<i>B. subtilis</i>	>3.0	15 min
EN1276	<i>P. aeruginosa</i>	>5.3	5 min	EN13697	<i>B. subtilis</i>	>3.0	15 min
EN1276	<i>E. coli</i>	>5.3	5 min	EN16615	<i>A. brasiliensis</i>	4.29 test area 1	2 min
EN1650	<i>C. albicans</i>	>4.3	15 min	EN16615	<i>B. subtilis</i>	3.18 test area 1	2 min

Product Application

Ready to use so can be applied directly to the surface ensuring even coverage. Physically remove contaminants from the surface using an appropriate cleanroom wipe and a recommended wiping technique for optimum contamination control. Leave for required contact time and wipe to dry. Contec HydroPure is very low residue, eliminating the need for a residue removal stage with alcohol or water.

Part No.	Description	Packaging
SBC056HP	Contec <i>Sterile</i> HydroPure 0.5L Capped	24 x 0.5L
SBT16HPW	Contec <i>Sterile</i> HydroPure 1L Trigger Spray	6 x 1L
SBC16HP	Contec <i>Sterile</i> HydroPure 1L Capped	6 x 1L
SBC56HP	Contec <i>Sterile</i> HydroPure 5L Capped	2 x 5L
FBC056HP	Contec HydroPure 0.5L Capped	24 x 0.5L
FBT16HP	Contec HydroPure 1L Trigger Spray	6 x 1L
FBC16HP	Contec HydroPure 1L Capped	6 x 1L
FBC56HP	Contec HydroPure 5L Capped	2 x 5L

Use biocides safely. Always read the label and product information before use.

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