




## Empfohlene Reinraumklassen

ISO 5|6|7|8|9

GMP A/B|C|D



## Tuch SterilSorb #C2-99IR/25

pure<sup>11</sup>-Nr.: 1106223, Marke: Contec

### Eigenschaften

- Steril
- Marke: Contec
- Flächengewicht in g/m<sup>2</sup>: 68 g/m<sup>2</sup>
- Anzahl Lagen: 1
- Struktur: Vliesstoff
- Material: Polyester/Zellulose
- Lieferform in Verpackung: gelegt
- Verpackungsform: Flachpackung (FlowPack)
- Anzahl in kleinster Unterverpackung: 25
- Besonders reißfest
- Chemikalienbeständigkeit
- Gebrauch: Einweg
- Hydroverfestigt
- Polyester/PES: 45 %
- Sterilisation validiert
- Zellulose: 55 %
- Flächengewicht in g/m<sup>2</sup>: 68 g/m<sup>2</sup>

## Material

- Polyester/Zellulose

## Verpackung

- 1500STK

## Produktvarianten

---

**pure<sup>11</sup>-Nr.: 11062232323WH, Tuch SterilSorb #C2-99IR/25**

Steril; Größe: 23 x 23 cm (9 x 9"); Farbe: Weiß / VE: 1500STK

---

**pure<sup>11</sup>-Nr.: 11062233131WH, Tuch SterilSorb #C2-1212IR**

Steril; Größe: 31 x 31 cm (12 x 12"); Farbe: Weiß / VE: 1050STK

---

# SterileSorb™ Wipes

## Cellulose/polyester nonwoven wipes

SterileSorb™ sterile hydroentangled polyester/cellulose nonwoven wipes are ideal for general purpose cleaning in sterile environments. Packaged in our sterile peelable outer bag and linear tear inner bag specifically designed for sterile products. Embossed lot numbers on both inner and outer bags ensure lot traceability. Sterile validation reports are available upon request.

SterileSorb is ideal for general purpose cleaning in biotech, pharmaceutical, and medical device manufacturing environments.




| Features                                                            | Benefits                                                                                                                                                                    |
|---------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Strong, durable hydroentangled polyester/cellulose wipes            | <ul style="list-style-type: none"> <li>• Very low levels of particles and fibers</li> <li>• Good abrasion and chemical resistance</li> <li>• Extremely absorbent</li> </ul> |
| Lot numbers on inner and outer bag                                  | <ul style="list-style-type: none"> <li>• Ensures lot traceability</li> </ul>                                                                                                |
| Validated sterile to a 10 <sup>-6</sup> SAL per ANSI/AAMI/ISO 11137 | <ul style="list-style-type: none"> <li>• Suitable for use in Grade A and B cleanrooms</li> </ul>                                                                            |
| Packaged in sterile peelable outer bags and linear tear inner bags  | <ul style="list-style-type: none"> <li>• Each bag is easy to open even when wearing gloves</li> <li>• Facilitates cleaning products entering isolators and RABS</li> </ul>  |

| Part No.   | Description                        |                                                                                     | Size                        | Packaging               |
|------------|------------------------------------|-------------------------------------------------------------------------------------|-----------------------------|-------------------------|
| C2-99IR/25 | SterileSorb Wipes,<br>Flat stacked |  | 9" x 9"<br>(230 x 230 mm)   | 25/bag;<br>60 bags/case |
| C2-1212IR  | SterileSorb Wipes,<br>Flat stacked |  | 12" x 12"<br>(305 x 305 mm) | 25/bag;<br>42 bags/case |

## Product Information

|                     |                                                                                                                                                                                                                                                                          |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Material            | • Cellulose/polyester                                                                                                                                                                                                                                                    |
| Construction        | • Hydroentangled                                                                                                                                                                                                                                                         |
| Packaging materials | • Outer bags (OB1, OB2), low density polyethylene (LDPE) <br>Case (CS), corrugated fiberboard (PAP)  |
| Environment         | • ISO 5 Grade A/B                                                                                                                                                                                                                                                        |



## Technical Data

| Attribute (units)                                              | Typical Value | Test Method                 |
|----------------------------------------------------------------|---------------|-----------------------------|
| Basis weight, nominal; (g/m <sup>2</sup> )                     | 68            | Contec Method               |
| Sorbent capacity; (mL/m <sup>2</sup> )                         | 314           | IEST-RP-CC004.3, Sec. 8.1   |
| Sorptive rate; (seconds)                                       | <1            |                             |
| Non-volatile residue, NVR                                      |               | IEST-RP-CC004.3, Sec. 7.1.2 |
| In deionized water; (g/m <sup>2</sup> )                        | 0.009         |                             |
| In isopropyl alcohol; (g/m <sup>2</sup> )                      | 0.009         |                             |
| Specific ions                                                  |               | IEST-RP-CC004.3, Sec. 7.2.2 |
| Sodium; (ppm)                                                  | 13.1          |                             |
| Chloride; (ppm)                                                | 9.5           |                             |
| Particles, readily releasable                                  |               |                             |
| Particles $\geq 0.5\mu\text{m}$ ; ( $\times 10^6/\text{m}^2$ ) | 32.7          | IEST-RP-CC004.2, Sec. 5.1   |
| Fibers $\geq 100\mu\text{m}$ ; ( $\times 10^3/\text{m}^2$ )    | 22.0          | IEST-RP-CC004.2, Sec. 5.2   |

## Recycle Symbols

PET 

HDPE 

LDPE 

PP 

PAP 

| Packaging  | EA/OB1 | OB1/OB2 | OB2/CS | EA/CS |
|------------|--------|---------|--------|-------|
| C2-99IR/25 | 25     | 6       | 10     | 1,500 |
| C2-1212IR  | 25     | 6       | 7      | 1,050 |

EA = each; OB = outer bag; CS = case

## Notes

- a) The data shown are typical values and should not be used as product specifications.  
b) Valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.  
c) Current and/or comparison data may be available. Please contact a Contec sales representative for details.

## Test Methods

1. CTM = Contec Test Method  
2. IEST-RP-CC004.3 = Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments, Institute of environmental Sciences and Technology, Rolling Meadows IL