

Distributor pure11



## Armschoner WeePro Steril #WL-MMS-00

pure<sup>11</sup>-Nr.: 1103514, Marke: Distributor pure11

### Eigenschaften

- Marke: Distributor pure11
- Gebrauch: Einweg
- Länge in cm: 50 cm
- Abschluss Unterarm: Gummiband und Daumenschlaufe
- Steril
- Anzahl in kleinster Unterverpackung: 2
- Farbe: Weiß

### Empfohlene

### Reinraumklassen

ISO 4|5|6|7|8|9

GMP A/B|C|D



## Material

- 

## Verpackung

- 100PAAR

## Produktvarianten

---

**pure<sup>11</sup>-Nr.: 1103514, Armschoner WeePro Steril #WL-MMS-00**

Steril; Farbe: Weiß; Länge: 50 cm / VE: 100PAAR

---

# MICROMIUM SERIES



**STERILE**

**WL-MMS-00**



## Key features

- New material **MICROMIUM®**, light and very low contaminant
- Cleanroom apparel compatible from Iso class 4
- Sealed packaging by pair with an Easy Tear system
- Overpack by 25 pairs
- Apparel sterilised by gamma irradiation

## Design features

- Bound seams and tunneled elastics
- Folded sleeve to gowning easier, and to prevent the contact with the outside surface
- Elastic cuffs and thumb loops
- Long length: 50cm

## Production setting

- Manufactured in a controlled environment and a facility holding ISO 9001:2008 certification

## Applications

- Pharmaceutical and cosmetic industries
- Biotech and research lab
- High-tech industries
- Medical device

## Packaging and tracking

100 pairs  
per carton



Big bag sealed  
by clip



Secondary transparent  
sealed bag by 25 pairs



One pair per sealed bag  
with «Easy Tear» system



Tracking number from carton to the garment

**FEATURES:** Sleeve with elastic top and wrist

Thumb loop

Elastics entunneled

**FABRICS AND SEAMS:** Microporous film

Bound seams

### PHYSICAL PERFORMANCE

Fabric Physical Properties	Test Method	Résultats	
Abrasion resistance	EN530	Class 2	
Flex cracking	ISO 7854 B	Class 6	
Trapezoidal Tear	ISO 9073-4	Class 1	
Tensile strength	ISO 13934-1	Class 1	
Puncture resistance	EN863	Class 1	

### FABRIC RESISTANCE TO PENETRATION BY LIQUIDS EN ISO 6530

Fabric Chemical Properties		Penetration	Repellency
Sulphuric Acid 30%		Class 3	Class 3
Sodium Hydroxide 10%		Class 3	Class 3

### STERILIZATION METHOD

GAMMA irradiation (25 kGy)

### CLEANLESS

HELMKE DRUM	Test method	Particule emission rate $\geq 0.3 \mu\text{M}$	Result
	IEST-RP-CC003.3	$\leq 2\ 000$	Category 1

### PACKAGING

TRIPLE PACKAGING	First	Second	Third
	2 units - Sealed and vacuum Easytear system	50 units closed by clip	200 units closed by clip

### SIZE

Length	Width		WeeSafe code
50 CM	20 CM		WL-MMS-00

### PACKING CARTON

	Units	Weight	Sizes	Volume
Carton	200	5.4 KG	43 X 34 X 35 CM	0,052 <sup>M3</sup>
Pallet	4 800	140 KG	80 X 120 CM	N/A

**GENCOD : 3701012802442**

**CUSTOM CODE : 6210109800**

WeeSafe – IMMEUBLE "LE FUJY"- 1 RUE DU CENTRE 93160 NOISY-LE-GRAND - FRANCE – Tél : + 33 (0)1 43 09 98 04  
Mail : [contact@weesafe.fr](mailto:contact@weesafe.fr)