




Tuch SterileWipe LP

pure¹¹-Nr.: 06211, Hersteller: ITW Texwipe



Zusammenfassung

- Neue pure11-Artikelnummer (ab 01.07.2023): 1106211
- Material: Polyester
- Gewicht: 155 g/m²
- Tränkung: Keine
- Gestrick aus 100% Polyester
- Vorgewaschen
- Hohe Absorptionskraft
- Gammasterilisiert (inkl. Sterilisationszertifikat)
- Lieferform: gelegt

Empfohlene Reinraumklassen

ISO ☐ 3 ☐ 4 ☒ 5 ☒ 6 ☒ 7 ☒ 8 ☒ 9

GMP ☐ ☐ ☒ A/B ☐ ☒ C ☒ D ☐

Produktvarianten

pure¹¹-Nr.: 062112323

Farbe: Weiß / Maße: 23 x 23 cm (9 x 9") / Herst.-Nr.: TX3211 / VE: 100 Stück

Quelle: <https://www.pure11.de/tuch-sterilewipe-lp>

SterileWipe™ LP

Knitted polyester wiper for cleaning critical sterile environments and dispensing solutions



Description

SterileWipe™ LP is a 100% continuous-filament polyester wiper that is low in particle generation and extractable levels. Each package is gamma irradiated to 10⁻⁶ Sterility Assurance Level.

SterileWipe™ LP provides excellent cleaning performance in critical sterile environments such as pharmaceutical aseptic fill areas, biotechnology manufacturing facilities, microbiological laboratories, sterile suites and prep rooms.

Features

- Gamma irradiated to 10⁻⁶ Sterility Assurance Level, according to AAMI Guidelines
- Lot specified information on each shipment simplifies record keeping:
 - Certificate of Processing confirming radiation dosage
 - Certification of Compliance
 - Expiration date marked on each package
- Cleanroom-laundered 100% polyester knit
- Independent quality control audits for sterility assurance

Benefits

- Withstands exposure to common cleaning and disinfecting solutions including bleaches, phenols and quaternary ammonium compounds
- Packaging permits alcohol wipedown of exterior bag before introduction of the inner bag into sterile suites
- Superior strength allows for rigorous wiping with minimal release of fibers and particles
- Wiper glides smoothly over surfaces
- Assures quality of individual lots and lot-to-lot traceability

Applications

- Developed to clean sterile environments without compromising the sterile nature of those environments
- Designed for cleaning and polishing equipment and environmental surfaces during and following the production flow
- Ideal for applying cleaning, sterilizing and disinfecting solutions

Products

TX Number	Description	Packaging
TX3211	SterileWipe™ LP 9" x 9" nominal (23 cm x 23 cm) polyester knit wipers	100 wipers/bag 10 inner bags of 10 wipers; 5 bags/case



North America
1210 South Park Drive
Kernersville, NC 27284
Tel (800) TEXWIPE
(336) 996-7046
Fax (336) 996-2297
www.texwipe.com
info@texwipe.com

Europe/Middle East
Skejby Nordlandsvej 307
DK-8200 Aarhus N
Denmark
Tel +45 87 400 220
Fax +45 87 400 222

Asia/Pacific
50 Tagore Lane
#02-01 Entrepreneur Centre
Singapore 787494
Tel +65 6468 9433
Fax +65 6468 6772

SterileWipe™ LP

TX3211

Performance Characteristics

Property	Typical Value	Test Method*
Basis weight	155 g/m ²	1, TM2
Absorbency		
Sorptive capacity	475 mL/m ²	1, TM3
Sorptive rate	<1 second	1, TM3

Contamination Characteristics

Property	Typical Value	Test Method*
Particles and fibers		
Particles 0.5-5.0 µm	10 x 10 ⁶ particles/m ²	1, 2, TM15
5.0-100 µm	403,000 particles/m ²	1, 2, TM15
Fibers: >100 µm	2,200 fibers/m ²	1, 2, TM15
Nonvolatile residue		
IPA extractant	0.08 g/m ²	1, TM1
DIW extractant	0.02 g/m ²	1, TM1
Ions		
Sodium	0.34 ppm	1, TM18
Potassium	0.17 ppm	1, TM18
Chloride	0.20 ppm	1, TM18

Note: The data in this table represent typical analyses of these wipers at the time of publication. These are not specifications.
ITW Texwipe continually refines both its processes and its products.

*Test Methods

- 1 – “Evaluating Wiping Materials Used in Cleanroom and Other Controlled Environments,” IEST-RP-CC 004.3, Institute for Environmental Sciences and Technology, Rolling Meadows, IL 2004; www.iest.org.
 - 2 – “Standard Method for Size-Differentiated Counting of Particles and Fibers Released from Clean Room Wipers Using Optical and Scanning Electron Microscopy,” E2090-00, ASTM International, West Conshohocken, PA, 2000; www.astm.org.
- TM – Refers to ITW Texwipe Test Method — available upon request, contact ITW Texwipe Customer Service at www.texwipe.com for a copy.

ITW Texwipe is the only wiper company to be ISO 9001:2000, 13485:2003, 14001:2004 and OHSAS 18001:1999 registered.